


Technical Evaluation & Validation of Locally Manufactured Ventilator (LMV)

Ventilator Type: Conventional

Applicant Organization/Team Name:

Conducted by




Military Institute of Science and Technology (MIST)

Mirpur Cantonment, Dhaka-1216



Version: MIST_Ventilator_TEV_1.0

[Handwritten signatures]



Summary of Technical Evaluation & Validation

Clause Ser.	Events	Remarks
A	Ventilation	
B	Gas and Electricity Supply	
C	Infection Control	
D	Monitoring and Alarm	
E	Biological Safety	
F	Software Safety	




G	Additional Observations	
---	-------------------------	--

Overall Recommendation

--

Signatures of the Technical Evaluation and Validation Committee



<p>.....</p> <p>Maj Md. Ashrafuzzaman, PhD, EME</p> <p>Dated.....</p>	<p>.....</p> <p>Lt Col Muhammad Nazrul Islam, PhD, Sigs</p> <p>Dated.....</p>	<p>.....</p> <p>Assistant Professor Golam Mostafa</p> <p>Dated.....</p>
<p>.....</p> <p>Name</p> <p>Dated.....</p>	<p>.....</p> <p>Name</p> <p>Dated.....</p>	<p>.....</p> <p>Name</p> <p>Dated.....</p>



Glossary

Before we proceed to the testing queries, this paper work would like to present elaborated introduction of some important abbreviated keywords:

Serial No	Abbreviated Form	Elaboration
1	ARDS	Acute Respiratory Distress Syndrome: a life-threatening form of respiratory failure where the lungs become severely inflamed due to an infection or injury and can't provide the body's vital organs with enough oxygen.
2	SIMV-PC	Synchronized Intermittent Mandatory Ventilation – Pressure Controlled: a mode of ventilation where the patient is allowed to take spontaneous breaths, the machine will assist the patients breathing when a spontaneous breath is taken. If the patient does not make a pre-set number of breaths a minute (i.e. 10) the machine provides mechanical ventilation to provide the set number.
3	CMV	Continuous Mandatory Ventilation
4	PCV	Pressure Controlled Ventilation
5	VCV	Volume Controlled Ventilation
6	PRVC	Pressure Regulated Volume Controlled: A mode of ventilation where a set tidal volume is delivered to the patient while maintain the lowest pressure possible in the airway, to avoid trauma.
7	CPAP	Continuous Positive Airway Pressure a non-invasive ventilation mode that provides a constant steady pressure to keep the lungs expanded.
8	BIPAP	Bi-level Positive Airway Pressure: a non-invasive ventilation mode that provides different levels of pressure when the patient inhales and exhales.
9	IPPV	Intermittent Positive Pressure Ventilation: a mandatory invasive ventilation mode used to replace a patient's breathing when they cannot for themselves. Can be either volume controlled or pressure controlled. It does not synchronize any patient breathing efforts.
10	PEEP	Positive End-Expiratory Pressure: The lower pressure applied to the patient's airway to allow them to breathe out, but not too much.
11	EPAP	Expiratory Positive Airway Pressure: Similar to PEEP, pressure applied to the airway on patient expiration to prevent collapse of the airway.
12	HMEF	Heat and Moisture Exchange Filter: device fitted to the patient end of the breathing system, contains hydrophobic medium that absorbs heat and moisture from the patients exhaled breath and uses absorbed moisture to humidify inhaled gases. Can also filter bacteria and viruses, this will be used on all patients. WARNING can affect delivered pressure.
13	RF	Radio Frequency: Many medical devices are sensitive to RF interference. Care should be taken to ensure that this is kept to a minimum.
14	EM	Electro Magnetic Emissions: Many medical devices are sensitive to EM interference. Care should be taken to ensure that this is kept to a minimum.
15	FiO ₂	Fraction of inspired oxygen: concentration of oxygen in the gas mixture that the patient inhales
16	AGSS	Anaesthetic gas scavenging system: where anaesthetic agents have been included in the gas mixture, this system is used to collect and remove exhaled gas to avoid exposure to health care professionals.
17	LMV	Locally Manufactured Ventilator




1. Introduction

COVID-19 pandemic caused by SARS-CoV-2 has recently increased the need of ventilators in hospitals throughout the world. Sudden increase in the demand of ventilators has created a crisis in the international market and many countries are now taking initiatives to produce locally designed ventilators. The whole world including Bangladesh is going through a critical period. Many individuals, independent researchers, research facilities, software firms, educational institutes and industries in Bangladesh have developed different categories of ventilators to save the life of the patient during the transition period. In order to cope with the limited resources, it is imperative to determine whether these ventilators are suitable for human use at all. This protocol aims at classifying locally made different types of ventilators based on their technical capabilities, efficiency, vulnerability and potential for application to humans. It follows the guidelines of intensive care medical professionals and anaesthetists to set a minimum criteria and performance of a ventilator to be called either as a Locally Manufactured Ventilator Emergency (LMV (E)) and Locally Manufactured Ventilator Conventional (LMV(C)). It is for alternative and conventional ventilators designed in an emergency situation in order to save life of a patient requiring urgent or long-term ventilation for respiratory failure caused by SARS-CoV-2 or any accident. Devices not meeting the standards of the protocol are likely to provide no benefit to the patient in need of emergency or long-term ventilation and might lead to increased harm.

2. Objectives

The objectives of the report is as follows-

1. Functionality testing of a locally made ventilator.
2. Evaluation of safety of a locally made ventilator.
3. Evaluation of viability for clinical testing.



3. Technical Evaluation & Validation

Although there are a number of factors to be monitored for the practicability of a Conventional ventilator, for the pandemic of COVID-19 only the following (minimal) specifications will be considered for its viability. This is to confer therapeutic benefit (urgent ventilation at acute condition) for a COVID-19 affected patient requiring invasive ventilation due to respiratory failure



Furthermore in the testing questionnaires' three words- 'Must', 'Should', and 'Could' have special definition as-

Must : Defines the minimum viable product clinically acceptable by clinicians



Should : Defines highly desirable features of considerable benefit for therapeutic use. As time is of the essence if omitting one of these features significantly accelerates development and production it should be considered

Could : Features or options often found in respirators, but are of significantly lower priority in terms of the current need and should not be considered if they delay production and development or the provision of more important features



Clause No	Items for Checking	Meets the Minimum Requirements (Yes/No)	Remarks
A. Ventilation			
1	1.1	Must have at least 1, optionally 2 modes of ventilation	
	1.2	Must have the CMV mode as- a. (Ideally) Pressure Regulated Volume Control, or b. pressure controlled ventilation (PCV) or c. Minimally a volume controlled ventilation (VCV)?	
	1.3	"Ideally PRVC/PCV, an adaptive mode where the tidal volume is set and the lowest possible pressure is delivered to achieve this volume. PCV where the user has to provide the adaptive control to achieve tidal volume is only acceptable if the tidal volume delivered is clearly displayed and the user can set patient specific upper and lower tidal volume alarms to alert to the need to adjust the pressure."	
	1.4	For the VCV test, the user sets a tidal volume and respiratory rate. The tidal volume is delivered during the inspiratory period. Acceptable only if additional pressure limiting controls are available.	





Clause No	Items for Checking	Meets the Minimum Requirements (Yes/No)	Remarks
1.5	The ventilation system should have a spontaneous breathing pressure support mode for those patients breathing to some extent themselves, e.g. BIPAP or SIMV-PC. The user sets an inspiratory pressure and an expiratory pressure. The ventilator can sense when a patient starts to breathe in and apply the inspiratory pressure, then sense when the patient starts to breathe out and apply the expiratory pressure (this pressure is still positive but lower than the inspiratory pressure)		
2	If a pressure support mode is provided the LMV (C) must failsafe automatically onto mandatory ventilation if the patient stops breathing in this mode.		
3	Inspiratory airway pressure, the higher pressure setting that is applied to make the patient breathe in	-	
	3.1 Plateau pressure should be adjusted to achieve volume and must be limited to 35 cm H ₂ O by default. It is acceptable if an option to increase this to 70 cm H ₂ O in exceptional circumstances is provided. This must require a positive decision and action by the user		
	3.2 Peak pressure should be no more than 2 cm H ₂ O greater than plateau pressure.		
	3.3 If volume control ventilation is used, the user must be able to set inspiratory airway pressure limit in the range at least 15 – 40 cm H ₂ O in at least increments of 5 cm H ₂ O.		
	3.4 There must be a mechanical failsafe valve that opens at 80 cm H ₂ O.		
4	PEEP: The pressure maintained in the breathing system during expiration. a. LMV(C) must provide a range 5-20 cm H ₂ O adjustable in 5 cmH ₂ O increments. b. The patient breathing system must remain pressurized to at least the PEEP level setting at all times.		
5	Inspiratory: Expiratory ratio (I:E). The proportion of each breathing cycle that is spent breathing in compared to breathing out. a. LMV(C) must provide 1:2.0 (i.e. expiration lasts twice as long as inspiration) as the default setting. b. LMV(C) could provide adjustable I:E in the range 1:1 – 1:3.		
6	Respiratory Rate: The number of breathing cycles every minute. LMV(C) must provide range 10 – 30 breaths per minute in increments of 2 (only in mandatory mode) that can be set by the user.		
7	Tidal Volume (V _T) setting, if provided. The volume of gas flowing into the lungs during one inspiratory cycle a. Must have at least one setting of 400ml +/- 10 ml. b. Should have 350ml and 450 ml options. c. Could have a range 250 – 600 ml in steps of 50ml. d. Could have a range up to 800 ml.		




Clause No	Items for Checking	Meets the Minimum Requirements (Yes/No)	Remarks	
B. Gas and Electricity				
1	Incoming Gas Supply-	-		
	1.1	All gas connectors and hoses must be medical grade.		
	1.2	Oxygen supply system and connections should be robust.		
	1.3	Average oxygen consumption must be no more than 6 lpm. This may be allowed to increase as greater certainty is gained over oxygen supply.		
	1.4	An LMV(C) may include an oxygen concentrator as the source of oxygen. Note these will typically be limited to 10 lpm 96% oxygen.		
2	Electricity Supply-	-		
	2.1	If mains powered, LMV(C) must connect to 220 V mains via standard 3pin plug.		
	2.2	If mains electricity is required for functioning, LMV(C) must have a battery backup of at least 20 minutes in case of mains electricity failure.		
	2.3	Could utilize hot swappable batteries so that it can be run on battery supply for an extended period, e.g. 2 hours for within hospital transfer.		
	2.4	Must avoid harmful RF or EM emissions that could interfere with other critical care equipment.		
3	Gas supply to patient-	-		
	3.1	User must be able to control inspired oxygen proportion (FiO ₂). Must provide a (50% or 60%) and 100% options.		
	3.2	Should provide control variable between 30 – 100 % in 10% steps.		
	3.3	Patient breathing system must contain HME filter.		
	3.4	The patient breathing system should avoid detrimental rebreathing.		
4	All elements in the gas pathway must meet biological safety and low-pressure oxygen safety standards, especially to minimize risk of fire or contamination of the patient's airway.			





Clause No	Items for Checking	Meets the Minimum Requirements (Yes/No)	Remarks
C. Infection control			
1	All parts coming into contact with the patient's breath must be both disposable and changeable.		
2	All working components of the device must be contained within an impermeable casing.		
3	All external surfaces must be cleanable in the likely event that they get respiratory secretions or blood splatter on them. Cleaning would be by healthcare workers manually wiping using an approved surface wipe with disinfectant or cloths and approved surface cleaning liquid.		
4	There will be a separately sourced HMEF-bacterial-viral filter between the machine and patient which may impact on resistance within the system, which may need to be accounted for with some designs. The pressure being delivered to the patient is the specified pressure. If the filter has a resistance of, say 2 cmH ₂ O at 30 lpm, the ventilator needs to output 37 cmH ₂ O to achieve a set 35 cmH ₂ O at the patient. This will need further detailed consideration. Viral filtering filters may have much higher resistance that may be clinically relevant.		
5	Could include facility for hot water humidifier to be included in breathing system		
D. Monitoring and Alarms			
Alarms, alarm limits, and priorities are complex areas to optimize for human usability. The key is to get enough alarms but not too many and for alarms to be clearly ranked so that more urgent patient safety problems are highlighted more. Early attention to this area is important, and should be built in from the start. To check the following points need some attentions-			
1	Must alarm at:	-	
	1.1 Gas or electricity supply failure.		
	1.2 Machine switched off while in mandatory ventilation mode.		
	1.3 Inspiratory airway pressure exceeded		
	1.4 Inspiratory and PEEP pressure not achieved (equivalent to disconnection alarm).		
2	Monitoring displayed continuously so the user can verify.	-	
	2.1 Must show the current settings of tidal volume, frequency, PEEP, FiO ₂ , ventilation mode.		
	2.2 Must show the actual current airway pressure		
	2.3 Should show the achieved tidal volume, breathing rate, PEEP, and FiO ₂ .		
	2.4 If pressure support mode is provided there must be real time confirmation of each patient breath and an alarm if below acceptable range.		
	2.5 Could provide CO ₂ monitoring.		



Clause No	Items for Checking	Meets the Minimum Requirements (Yes/No)	Remarks
E. Biological Safety			
1	Materials of Construction (raw materials) a. The materials of construction must be medical grade and biocompatible. b. Polyvinyl chloride (PVC) must be avoided in the patient gas pathway c. PVC should be avoided elsewhere		
2	Hazard Mitigation:	-	
	2.1 Leachable substances (in condensate): chemical removed from the medical device by the action of water, other liquids or other gases related to the use of the medical device- insure a HME filter is used between the ventilator and breathing system		
	2.2 If any glue or adhesive is used, it must be medical grade and should not emit Volatile Organic Compounds (VOC) or leach chemicals during operation.		
F. Software Safety			
	Software in a high-risk device like LMV(C) will almost certainly have the capability to cause serious injury or death if risk control measures are not adequately implemented.	-	
1	Are the features claimed in the software requirement specification (SRS) implemented into the software developed?		
2	Are the ways of exception handling clearly addressed in the software in the software architecture and software design?		
3	Is the verification and validation process planned, conducted (in-house testing) and reported?		
4	Is the user interface of the software user-friendly? a. Visual Alarm Signals as well as audible alarm signals b. Labels clearly visible from at least 4ft distance c. Access distance to set up the RR, TV and I:E ratio is not more than 2 clicks.		
G. Miscellaneous			
1	Must be reliable. LMV(C) must be capable of continuous operation (100% duty cycle) for 12 Hours.		
2	Could be floor standing and light enough to mount on patient bed with orientation independent functioning.		
3	Should be as robust as possible.		





Clause No	Items for Checking	Meets the Minimum Requirements (Yes/No)	Remarks
4	It must be intuitive to use for qualified medical personnel, but these may not be specialists in ventilator use. a. Must include Instructions for Use. b. Instructions for use should be built into the labeling of the ventilator, e.g. with 'connect this to wall' etc. c. Must include clear labeling of all critical functions and controls using standard terms, pictograms and colors that will be readily recognized by healthcare staff.		
5	Must Provide User manual		
6	Service Manual should provide and is highly recommended. For end user delivery it is a must.		
H. Additional Observations			
1			
2			
3			
4			
5			
6			
7			
8			
9			
10			




4. Ventilation Test

4.1 Volume Controlled Ventilation Test (Compliance)

Test No	Test Lung Compliance ml/cmH ₂ O ± 10%	Test Lung Resistance cmH ₂ O/l/s ± 10%	Tidal Volume ml	Rate Min ⁻¹	I:E	O ₂ %	PEEP	Findings	Remarks
1	50	5	500	20	1:2	50-60	5		
2	50	5	500	20	1:2	90-100	5		
3	50	5	500	12	1:2	50-60	5		
4	50	5	500	12	1:2	90-100	5		
5	50	5	500	20	1:2	50-60	10		
6	50	5	500	20	1:2	90-100	10		
7	50	5	500	12	1:2	50-60	10		
8	50	5	500	12	1:2	90-100	10		
9	50	5	500	20	1:2	50-60	15		
10	50	5	500	20	1:2	90-100	15		
11	50	5	500	12	1:2	50-60	15		
12	50	5	500	12	1:2	90-100	15		
13	20	5	500	20	1:2	50-60	5		
14	20	5	500	20	1:2	90-100	5		
15	20	5	500	12	1:2	50-60	5		
16	20	5	500	12	1:2	90-100	5		
17	20	5	500	20	1:2	50-60	10		
18	20	5	500	20	1:2	90-100	10		
19	20	5	500	12	1:2	50-60	10		
20	20	5	500	20	1:2	50-60	10		
21	20	5	500	20	1:2	50-60	15		
22	20	5	500	20	1:2	90-100	15		
23	20	5	500	12	1:2	50-60	15		
24	20	5	500	20	1:2	50-60	15		





25	10	5	500	20	1:2	50-60	5		
26	10	5	500	20	1:2	90-100	5		
27	10	5	500	12	1:2	50-60	5		
28	10	5	500	12	1:2	90-100	5		
29	10	5	500	20	1:2	50-60	10		
30	10	5	500	20	1:2	90-100	10		
31	10	5	500	12	1:2	50-60	10		
32	10	5	500	20	1:2	50-60	10		
33	10	5	500	20	1:2	50-60	15		
34	10	5	500	20	1:2	90-100	15		
35	10	5	500	12	1:2	50-60	15		
36	10	5	500	20	1:2	50-60	15		





4.2 Volume Controlled Ventilation Test (Resistance)

Test No	Test Lung Compliance ml/cm H ₂ O ± 10%	Test Lung Resistance cm H ₂ O/l/s ± 10%	Tidal Volume ml	Rate Min ⁻¹	I:E	O ₂ %	PEEP	Findings	Remarks
1	50	5	500	20	1:2	50-60	5		
2	50	5	500	20	1:2	90-100	5		
3	50	5	500	12	1:2	50-60	5		
4	50	5	500	12	1:2	90-100	5		
5	50	5	500	20	1:2	50-60	10		
6	50	5	500	20	1:2	90-100	10		
7	50	5	500	12	1:2	50-60	10		
8	50	5	500	12	1:2	90-100	10		
9	50	5	500	20	1:2	50-60	15		
10	50	5	500	20	1:2	90-100	15		
11	50	5	500	12	1:2	50-60	15		
12	50	5	500	12	1:2	90-100	15		
13	20	20	500	20	1:2	50-60	5		
14	20	20	500	20	1:2	90-100	5		
15	20	20	500	12	1:2	50-60	5		
16	20	20	500	12	1:2	90-100	5		
17	20	20	500	20	1:2	50-60	10		
18	20	20	500	20	1:2	90-100	10		
19	20	20	500	12	1:2	50-60	10		
20	20	20	500	20	1:2	50-60	10		
21	20	20	500	20	1:2	50-60	15		
22	20	20	500	20	1:2	90-100	15		
23	20	20	500	12	1:2	50-60	15		
24	20	20	500	20	1:2	50-60	15		
25	10	50	500	20	1:2	50-60	5		
26	10	50	500	20	1:2	90-100	5		

Signature





27	10	50	500	12	1:2	50-60	5		
28	10	50	500	12	1:2	90-100	5		
29	10	50	500	20	1:2	50-60	10		
30	10	50	500	20	1:2	90-100	10		
31	10	50	500	12	1:2	50-60	10		
32	10	50	500	20	1:2	50-60	10		
33	10	50	500	20	1:2	50-60	15		
34	10	50	500	20	1:2	90-100	15		
35	10	50	500	12	1:2	50-60	15		
36	10	50	500	20	1:2	50-60	15		





4.3 Volume Controlled Ventilation Test (Tidal Volume)

Test No	Test Lung Compliance ml/cm H ₂ O ± 10%	Test Lung Resistance cm H ₂ O/l/s ± 10%	Tidal Volume ml	Rate Min ⁻¹	I:E	O ₂ %	PEEP	Findings	Remarks
1	50	5	300	20	1:2	50-60	5		
2	50	5	300	20	1:2	90-100	5		
3	50	5	300	12	1:2	50-60	5		
4	50	5	300	12	1:2	90-100	5		
5	50	5	300	20	1:2	50-60	10		
6	50	5	300	20	1:2	90-100	10		
7	50	5	300	12	1:2	50-60	10		
8	50	5	300	12	1:2	90-100	10		
9	50	5	300	20	1:2	50-60	15		
10	50	5	300	20	1:2	90-100	15		
11	50	5	300	12	1:2	50-60	15		
12	50	5	300	12	1:2	90-100	15		
13	20	20	300	20	1:2	50-60	5		
14	20	20	300	20	1:2	90-100	5		
15	20	20	300	12	1:2	50-60	5		
16	20	20	300	12	1:2	90-100	5		
17	20	20	300	20	1:2	50-60	10		
18	20	20	300	20	1:2	90-100	10		
19	20	20	300	12	1:2	50-60	10		
20	20	20	300	20	1:2	50-60	10		
21	20	20	300	20	1:2	50-60	15		
22	20	20	300	20	1:2	90-100	15		
23	20	20	300	12	1:2	50-60	15		
24	20	20	300	20	1:2	50-60	15		
25	10	50	300	20	1:2	50-60	5		
26	10	50	300	20	1:2	90-100	5		

Amman





27	10	50	300	12	1:2	50-60	5		
28	10	50	300	12	1:2	90-100	5		
29	10	50	300	20	1:2	50-60	10		
30	10	50	300	20	1:2	90-100	10		
31	10	50	300	12	1:2	50-60	10		
32	10	50	300	20	1:2	50-60	10		
33	10	50	300	20	1:2	50-60	15		
34	10	50	300	20	1:2	90-100	15		
35	10	50	300	12	1:2	50-60	15		
36	10	50	300	20	1:2	50-60	15		





4.4 Pressure Controlled Ventilation Test (15 cm H₂O)

Test No	Test Lung Compliance ml/cm H ₂ O ± 10%	Test Lung Resistance cm H ₂ O/l/s ± 10%	Plateau Pressure cm H ₂ O	Rate Min ⁻¹	I:E	O ₂ %	PEEP	Findings	Remarks
1	50	5	15	20	1:2	50-60	5		
2	50	5	15	20	1:2	90-100	5		
3	50	5	15	12	1:2	50-60	5		
4	50	5	15	12	1:2	90-100	5		
5	50	5	15	20	1:2	50-60	10		
6	50	5	15	20	1:2	90-100	10		
7	50	5	15	12	1:2	50-60	10		
8	50	5	15	12	1:2	90-100	10		
9	50	5	15	20	1:2	50-60	15		
10	50	5	15	20	1:2	90-100	15		
11	50	5	15	12	1:2	50-60	15		
12	50	5	15	12	1:2	90-100	15		
13	20	20	15	20	1:2	50-60	5		
14	20	20	15	20	1:2	90-100	5		
15	20	20	15	12	1:2	50-60	5		
16	20	20	15	12	1:2	90-100	5		
17	20	20	15	20	1:2	50-60	10		
18	20	20	15	20	1:2	90-100	10		
19	20	20	15	12	1:2	50-60	10		
20	20	20	15	20	1:2	50-60	10		
21	20	20	15	20	1:2	50-60	15		
22	20	20	15	20	1:2	90-100	15		
23	20	20	15	12	1:2	50-60	15		
24	20	20	15	20	1:2	50-60	15		
25	10	50	15	20	1:2	50-60	5		
26	10	50	15	20	1:2	90-100	5		
27	10	50	15	12	1:2	50-60	5		





28	10	50	15	12	1:2	90-100	5		
29	10	50	15	20	1:2	50-60	10		
30	10	50	15	20	1:2	90-100	10		
31	10	50	15	12	1:2	50-60	10		
32	10	50	15	20	1:2	50-60	10		
33	10	50	15	20	1:2	50-60	15		
34	10	50	15	20	1:2	90-100	15		
35	10	50	15	12	1:2	50-60	15		
36	10	50	15	20	1:2	50-60	15		



4.5 Pressure Controlled Ventilation Test (30 cm H₂O)

Test No	Test Lung Compliance ml/cm H ₂ O ± 10%	Test Lung Resistance cm H ₂ O/l/s ± 10%	Plateau Pressure cm H ₂ O	Rate Min ⁻¹	I:E	O ₂ %	PEEP	Findings	Remarks
1	50	5	30	20	1:2	50-60	5		
2	50	5	30	20	1:2	90-100	5		
3	50	5	30	12	1:2	50-60	5		
4	50	5	30	12	1:2	90-100	5		
5	50	5	30	20	1:2	50-60	10		
6	50	5	30	20	1:2	90-100	10		
7	50	5	30	12	1:2	50-60	10		
8	50	5	30	12	1:2	90-100	10		
9	50	5	30	20	1:2	50-60	15		
10	50	5	30	20	1:2	90-100	15		
11	50	5	30	12	1:2	50-60	15		
12	50	5	30	12	1:2	90-100	15		
13	20	20	30	20	1:2	50-60	5		
14	20	20	30	20	1:2	90-100	5		
15	20	20	30	12	1:2	50-60	5		
16	20	20	30	12	1:2	90-100	5		
17	20	20	30	20	1:2	50-60	10		
18	20	20	30	20	1:2	90-100	10		
19	20	20	30	12	1:2	50-60	10		
20	20	20	30	20	1:2	50-60	10		
21	20	20	30	20	1:2	50-60	15		
22	20	20	30	20	1:2	90-100	15		
23	20	20	30	12	1:2	50-60	15		
24	20	20	30	20	1:2	50-60	15		
25	10	50	30	20	1:2	50-60	5		
26	10	50	30	20	1:2	90-100	5		



27	10	50	30	12	1:2	50-60	5		
28	10	50	30	12	1:2	90-100	5		
29	10	50	30	20	1:2	50-60	10		
30	10	50	30	20	1:2	90-100	10		
31	10	50	30	12	1:2	50-60	10		
32	10	50	30	20	1:2	50-60	10		
33	10	50	30	20	1:2	50-60	15		
34	10	50	30	20	1:2	90-100	15		
35	10	50	30	12	1:2	50-60	15		
36	10	50	30	20	1:2	50-60	15		

Documents Attached with this Report

- 1.
- 2.
- 3.
- 4.
- 5.

.....END OF REPORT.....

[Handwritten signatures and initials]